

# Dental high-speed handpiece

## LOTUS 401(P)

ORDER CODE:T-009

Push Button Type



### DESCRIPTION

Standard head high-speed handpiece  
Push button type  
Single spray system

### SUPPLY AIR PRESSURE

177-216Kpa(1.8-2.2kgf/cm<sup>2</sup>)

### ROTATION SPEED

380, 000-420,000r/min

### CHUCK

Push button type

### NOISE

65-69dB

### Compatible

Lotus 401(P)cartridge can be used for Lotus 401 handpieces and NSK PAN AIR 2000 handpieces.

### Upgradable

NCK PANA AIR (PA-S) handpiece&Being's primary Lotus 401 handpiece can be upgraded to push button type by changing the cartridge and the head cap into that of Lotus 401(P)

### Comfortable handle

Slide-proof, comfortable to hold and easy to clean

### Push button type

It is convenient to change the bur with the push button

#### ★ RELEVANT SPARE PARTS:



401M back cap



401M cartridge



401M cap spanner