

# DENTAL LOW SPEED HANDPIECE

## ROSE 201-M4

**ORDER CODE:T-010**

Specially for dental clinical treatment and prevention.

ROSE 201-CA  
**ORDER CODE:T-011**



- E-type connector
- Precise technology
- Exquisite appearance
- Durable and autoclavable
- Easy to assemble

ROSE 201-SH  
**ORDER CODE:T-012**



ROSE 201AM-M4  
**ORDER CODE:T-013**



ROSE 201AM-B2  
**ORDER CODE:T-014**



## Specifications

Air-in Pressure(kg/cm <sup>2</sup> )	Motor Rotation(rpm)	Air Wastage:(N1/min)
2.5	22,000	42
3	25,000	51
4	27,000	72

\*Speed may vary slightly depending on the back-end configuration, the direction of rotation, and type of hose used.

## Relevant Spare Parts:

